


Establish a  
nephrostomy tract  
with precision  
and ease.



**Ultraxx™**  
NEPHROSTOMY BALLOON



[www.cookmedical.com](http://www.cookmedical.com)

# Ultraxx™

**NEPHROSTOMY BALLOON**

The Ultraxx™ Nephrostomy Balloon is used to dilate the musculofascia, establishing a nephrostomy tract during a percutaneous procedure. The balloon has a rated burst pressure of 20 atm, and the functional dilating portion of the balloon is 15 cm long. The combination of high pressure and balloon length helps a physician establish a nephrostomy tract in the most demanding cases.

The Ultraxx is available with either a standard PTFE or a clear vinyl access sheath.



The Ultraxx is available for purchase as a set with the Cook Inflation Device.

## Features

---

Port grip                      Provides added control and aids handling during the inflation process.

---

Sheath lock                    Holds the sheath in place and out of the way during dilation.

---

Dual inflation/  
deflation ports                Provide quick, uniform pressure during inflation. When the balloon is deflated, the ports provide complete fluid evacuation for an uncomplicated removal process.

---

Platinum band                Provides maximum radiopacity to ensure precise balloon dilation positioning in even the most demanding cases.


---

Tapered Tip                    Eases insertions, requires minimal predilation and allows placement over a .038 inch wire guide.

---

Clear sheath                    Provides improved visualization.

---



## Ultraxx™

NEPHROSTOMY BALLOON

Ultraxx features a high pressure balloon for establishing and maintaining a nephrostomy tract.

Order Number	Reference Part Number	Inflated Balloon Diameter mm	Balloon Length cm
<b>Set</b>			
With Clear Vinyl Sheath			
G46037	UNBS-6-15-CS	6	15
G46038	UNBS-8-15-CS	8	15
G30357	UNBS-10-15-CS	10	15
With PTFE Sheath			
G46036	UNBS-6-15	6	15
G46024	UNBS-8-15	8	15
G30365	UNBS-10-15	10	15

Order Number	Reference Part Number	Inflated Balloon Diameter mm	Balloon Length cm
<b>Balloon Catheter Only</b>			
With Clear Vinyl Sheath			
G46023	UNB-6-15-CS	6	15
G46035	UNB-8-15-CS	8	15
G30358	UNB-10-15-CS	10	15
With PTFE Sheath			
G46033	UNB-6-15	6	15
G46034	UNB-8-15	8	15
G30364	UNB-10-15	10	15

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

### Customer Service Centers - Urology

**EMEA:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
 Poland: +48 223060159, [pl.orders@cookmedical.com](mailto:pl.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

[www.cookmedical.com](http://www.cookmedical.com)

**Americas:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
 E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)  
**Australia:**  
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283  
 E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)

